


MRI Simulation

MSc Martha Rebeca Canales Fiscal.

RF effect

- Open the simulator: <https://www.drcmr.dk/BlochSimulator/>
- You start by looking at the magnetization of a sample represented with a gray bar. To begin with the activity, be sure the specifications are as follows:



Mode: Equilibrium

Relaxation: off

B0: 2


B1: 0

B1freq: not relevant at this point (it can have any value)


Gradients: both 0

Frame: Stationary.


- You are going to look at the effect of applying an RF (B_1 field).
- Click on View and be sure all the options are activated.



- Click on Field and be sure B_0 is on 3T, let B_1 freq in any value. Then put B_1 in 0.3 and look at the effect on the magnetization (grey bar).




- You can see the magnetization precessing arbitrarily.
- To see more clearly the effect you should change the frame of reference. Observe how much and in which manner the magnetization moved away from the original direction.
- Click on Frame and choose the frame of reference of the field B0.




- The effect is not too strong because the frequency of B1 is not the resonant frequency. **Change the value of B1freq until the torque points towards only one direction. Once you find the appropriate frequency register it for your report**

The signal that we measure from the precessing magnetization is proportional to the size of the transversal magnetization (R_F/B_1). It can be seen in the length of the shadow (marked in red).




- Change between B0 and stationary frames of reference to observe by yourself and understand how the complete phenomenon is happening.
- Once you have a strong signal (the length of the shadow in large) turn off B1.




Relaxation times

- Until now, you explored the effect that the transversal field has on the magnetization. Now you should return it to the equilibrium. For doing so, you need to change the relaxation times.
- First look at the effect of the T2 relaxation time.
- Click on Relaxation and change T2 to 5 seconds.



Observe how the
transversal
magnetization decays

- Now change T1 to 5 seconds.
- Observe how it returns to the equilibrium state.




Complete process

- Now repeat the entire process by yourself using different relaxation times combinations.
- 1) Start in the equilibrium.
 - 2) Turn on B1.
 - 3) Set the resonant frequency.
 - 4) When you obtain a strong signal, turn off B1.
 - 5) Use different relaxation times and observe the effect when changing the frames of reference from Stationary to B0 and again to Stationary as many times as you need to understand the process (at least once).

Remember to take ss during convenient times for your report.


Different tissues

- Change the mode to Mixed matter. You are going to see the magnetizations corresponding to samples of different tissues with different colors. At the beginning all of them are together so you can't differentiate them.




The magnetizations of all the tissues are here all together.

- Now repeat the process that you previously learned and observe how the magnetizations of the different tissues respond differently (**Take evidence with ss**). **Explain in your report the differences (signal's strength and time) between the magnetizations referring to them by colors.**



Spin Echo


- Now you are going to explore spin-echo and with it the dephasing effect.
- Start with the specifications of the image. Note that Inhomogeneity mode was chosen. The commands in the bottom allow you to apply a quick pulse of 90° in the x-direction and 180° in the y-direction.





- Apply a 90°_x hard pulse.



- As the sample is not homogeneous you can see how the magnetizations are dephasing. **Take an ss of this effect.**
- After the magnetizations are dephased and they start to return to the equilibrium state apply a 180°_y pulse.



- Note how you obtain an echo with this process. If you do it at the correct timing the echo is going to be strong. **Try to obtain your best result and take a ss of the graph.** Play with different specifications (B_0 , T_1 , T_2 , etc.).



MRI 2nd Activity

Team SPECT

Mei Li Luisa Cham Perez


Ana Lucía Soria Cardona

Graciela Alejandra Rincón López

Natalia Verónica Flores Del Río


Marcela Enriquez López

step 1



step 2

Selección de las view y se establecieron los rangos de B0, B1 y B1freq



step 3

▸ Bloch Simulator

▸ Relaxation: Off

▸ View: B1,Torque,Mx,|Mxy|,Mz

▾ Fields

B0		3
B1		0.3
B1freq		3

▸ Gradients: $G_x=0$, $G_y=0$


▸ Frame: B0



step 4

Frequency at which the torque only moves at one direction. B1Freq= 2.8

B1Freq= 2.8



Our example

$$BO = 4$$

$$B1 = 0.9$$

$$B1Freq = 4$$

Case

▸ Bloch Simulator

▾ Relaxation

T1 Infinity

T2 Infinity

▸ View: B1,Torque,Mx,|Mxy|,Mz

▾ Fields



B0 4

B1 0.9

B1freq 4

Gradients: $G_x=0$, $G_y=0$

Frame: B0



Equilibrium

90° hard

90° soft

180°

Spoil

Non-rep. exc.



[Share](#) [in](#) [whatsapp](#) [twitter](#) [facebook](#) [telegram](#)

▸ Bloch Simulator

▾ Relaxation

T1 Infinity

T2 Infinity

▸ View: B1,Torque,Mx,|Mxy|,Mz

▾ Fields



B0 4

B1 0

B1freq 4

▸ Gradients: $G_x=0$, $G_y=0$

▸ Frame: B0



Equilibrium

90° hard

90° soft

180°

Spoil

Non-rep. exc.



▸ Bloch Simulator

▾ Relaxation

T1	<input type="text" value="9"/>
T2	<input type="text" value="9"/>



▸ View: B1,Torque,Mx,|Mxy|,Mz

▾ Fields

B0	<input type="text" value="4"/>
B1	<input type="text" value="0"/>
B1freq	<input type="text" value="4"/>

▸ Gradients: $G_x=0$, $G_y=0$

▸ Frame: B0



Lucía Soria

Equilibrium ▾

90°_x hard ▾

90°_x soft ▾


180°_y ▾

Spoil ▾

Non-rep. exc. ▾

||

Mixed matter




With the B1 frame, all the components are moving in the same direction.

Due to the different properties found in the simulation there are affected in different ways by the ultramagnetic field. This is seen as in the simulation the pieces move towards the same direction, yet they move at a different speed.


Looking at the different colors of pieces we can observe delays on the rotation of them, it can be seen how the green piece is moving at a faster rate than the blue piece.

Mixed matter

When editing the B1 and the relaxation, a faster decrease was noted by the white and blue bars, while the green one kept a relatively good size and rotation time.



90° hard angle



In the 90° hard angle the magnetizations are dephasing and we can appreciate that the echo is greater than m_x but is not strong enough to be greater than m_z .

180°_y angle

When the magnetizations are returning to the equilibrium state the 180°_y angle is applied and we can observe that the echo is stronger than both m_x and m_z .

